

P O W E R C H A I R



Seat Colors
Charcoal Vinyl

Body Colors
Candy Apple Red

Blue

Roll

Alanté Jr.
BY GOLDEN

Featuring a 250 lb. weight capacity, the Alanté Jr.'s new standard captain's seat offers a comfortable seating option on a compact, rear-wheel drive base. Fully transportable, the Alanté Jr. disassembles quickly and easily into three pieces.



Rise.

Roll.

Rest.

Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

Alanté GP-200 Specifications

| | |
|--|-----------------------|
| Model | GP200 |
| Medicare Code | K0823 |
| Weight Capacity | 250 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 3.5 mph |
| Operating Range ¹ | 9.2 miles |
| Ground Clearance to transaxle/motor | 1.5" |
| Turning Radius | 27" |
| Type Batteries | 2-22AH |
| Colors | Red, Blue |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Charger | Off Board |
| Controller | PG VR-2 40 Amp Remote |
| Dimensions: | |
| Length | 32.5" |
| Width | 23" |
| Height (ground to top of back of seat) | 35.5" - 36.5" |
| Height (ground to top of back of headrest) | 42" - 46" |
| Ground to Top of Seat | 20.5" - 21.5" |
| Number of Seat Height Positions | 2 |
| Size of Increments | 1" |
| Weight of Unit: (assembled) | |
| Base | 56 lbs |
| Seat w/Arms | 38 lbs |
| Batteries | 30 lbs/pack |
| Tires: | |
| Rear | 9" |
| Casters | 7" |
| Rear Anti-Tippers | 2.5" |
| Standard Seat: | |
| Back Height - no headrest | 18.5" |
| Back Height with headrest | 25" - 28" |
| Width x Depth | 18" X 16" |
| Color | Gray Vinyl |
| Footrest | |
| Height adjustable | No |
| Angle adjustable | Yes |
| Warranty: | |
| Frame | Lifetime |
| Drive Train | 1 Year |
| Electronics ² | 1 Year |

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



Disassembles quickly and easily into three pieces!



Unique trigger release for easy battery pack removal

Features

- Low center of gravity
- Excellent maneuverability
- 27" turning radius
- 250 lb weight capacity
- Gorgeous colors

Lifetime Warranty:
Power Chair Frame

One Year Limited Warranty:
Electronic Controller
Drive Train Components

Six-Month Limited Warranty:
Plastic Parts, Rubber Parts, Bearings
(Excluding plastic shroud)

Power Chair Accessories and Options

- Cane Holder
- Lap Belt
- Elevating Leg Rests