

P O W E R C H A I R



Seat Colors  
Charcoal/Gray Vinyl



Coffee/Sand Vinyl



Candy Apple Red



Onyx



Body Colors  
Blue  
Hunter Green



# Roll

golden  
**COMPASS**™

True, stable platform, center-wheel drive power chair with an articulating chassis for the ultimate in maneuverability and stability.



**GOLDEN**  
technologies

YOUR LIFE IN MOTION

## Golden Compass™ Specifications

Model Number	GP600CC	GP600SS
Medicare Code	K0823	K0822
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Center	Center
Maximum Speed	5 mph	5 mph
Operating Range <sup>1</sup>	20 miles	20 miles
Ground Clearance to transaxle/motor	3.5"	3.5"
Turning Radius	19"	19"
Type Batteries	2-U1 or 2-NF22	2-U1 or 2-NF22
Colors	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	Dynamic Shark 60 Amp Remote	Dynamic Shark 60 Amp Remote
<b>Dimensions<sup>3</sup>:</b>		
Length	38.5"	38.5"
Width	24.5"	24.5"
Height (ground to top of back of seat)	38.25" - 42.75"	38.25" - 42.75"
Height (ground to top of back of headrest)	45.25" - 53.75"	45.25" - 53.75"
Ground to Top of Seat	19.75" - 24.25"	17.75" - 22.25"
Number of Seat Height Positions	10	10
Size of Increments	0.5"	0.5"
<b>Weight of Unit: (assembled)</b>	189 lbs	186 lbs
Base	92 lbs	92 lbs
Seat w/Arms	49 lbs	46 lbs
Batteries	22 lbs ea.	22 lbs ea.
<b>Tires:</b>	Flat Free	Flat Free
Rear	10"	10"
Casters	6"	6"
<b>Standard Seat<sup>3</sup>:</b>	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 32.5"
Width x Depth	16" X 16", 18" X 18", 20" X 20"	16" X 16", 18" X 18", 20" X 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
<b>Optional Seats<sup>3</sup>:</b>	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	21"	23"
Back Height with headrest	28.5" - 32.5"	30.5" - 32.5"
Width	16", 18", 20", 22" or 24"	16", 18", 20", 22" or 24"
Depth	16", 18" or 20"	16", 18" or 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
<b>Optional Seats 2<sup>3</sup>:</b>	Stadium Style	Pro-Flex Seat
Back Height - no headrest	15.5"	N/A
Width	18"	16", 18", 20", or 22"
Depth	16"	16", 18" or 20"
Color	Black Vinyl	Black
<b>Seat Options:</b>		
Pneumatic Rise Seat Option	N/A	N/A
Electric Rise Seat Option	N/A	N/A
<b>Footrest</b>	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
<b>Warranty:</b>		
Frame	Lifetime	Lifetime
Drive Train	1 Year	1 Year
Electronics <sup>2</sup>	1 Year	1 Year



### 4 exciting colors



Candy Apple Red



Blue



Hunter Green



Onyx

### Power Chair Accessories and Options

- Stadium seat
- High back seat
- Pro Flex seat
- Van/pan seat
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Positioning belt
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go

<sup>1</sup> Battery range will vary due to rider weight, drive surface and terrain.

<sup>2</sup> Electronics warranty excludes batteries.

<sup>3</sup> There is a +1"/-1" tolerance on all dimensions.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.