

FULL SIZE PORTABLE SCOOTER



Seat Color
Black Vinyl



Body Colors
Blue
Red



Ride

GOLDEN
LITWAY™

4 wheel jockey style

The Liteway is in a class by itself! Experience the latest design in transport scooters with its rugged tubular frame, easy battery access and a unique storage basket under the seat! With plenty of legroom and a plush, comfortable seat, the Liteway is available in Red and Blue.



Rise.

Roll.

Rest.

Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

4 wheel jockey style Specifications

Model:	GL-100
Length:	41"
Width:	22"
Tires:	Foam-filled
Front Tire:	7"
Rear Tire:	9"
Weight Capacity:	300 lbs.
Max Speed:	4 mph
Ground Clearance:	3"
Battery Range using 2 15Ah batteries at Max Weight:	7.5 miles*
Turning Radius:	32"
Seat Weight:	21 lbs.
Battery Pack Weight:	19 lbs.
Front Half Weight:	39 lbs.
Rear Half Weight:	25 lbs.
Total weight including batteries:	106 lbs.
Freewheel Mode:	Yes
Electronic Speed Control:	Yes
Brakes:	Electro mechanical
Charger:	Off-board
Rear-wheel Drive:	Sealed transaxle motor, 24V DC
Controller:	45 amp non-programmable
Warranty – frame:	Limited lifetime
Warranty – electronics:	2 years
Warranty – drive train:	3 years



Liteway shown with optional seat.



The unique under-seat basket provides a large cargo area that's safe and secure.

Features

- 300 lb. weight capacity
- Rugged anti-rust steel frame tubing
- Rugged rubber floor mat for optimum durability
- Disassembles quickly and easily
- Dual throttle controls with both thumb and finger levers



Remove the seat and fold the tiller flat for easy storage in any trunk!



The Liteway comes standard with a plush, comfortable swivel seat and flip-up armrests for easy transfers.

Blue



Red



* Battery range will vary due to rider weight, drive surface and drive slope.

** Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.