

Evaluated By: _____

Company: _____

Mark For: _____

Address: _____

Date: _____

Account #: _____

Phone: _____

GC240 - COMPANION SCOOTER - (3) Three Wheel		K0806
<input checked="" type="checkbox"/>	GC240	COMPANION - 3 Wheel Scooter
	300 lb Capacity	\$2,735.00
with 18" X 18" Black Vinyl Stadium Seat		

Totals
\$

STANDARD COLORS - Please Choose One!

Red

 Blue

OPTIONAL SEATING - Additional Charges Apply

Please Choose Color	Stadium Style Seat Sizes																				
<input type="checkbox"/> Charcoal Vinyl	<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width: 10%;"></td> <td style="width: 10%; text-align: center;">Select One!</td> <td colspan="3" style="text-align: center;">Width</td> </tr> <tr> <td></td> <td></td> <td style="text-align: center;">16"</td> <td style="text-align: center;">18"</td> <td style="text-align: center;">20"</td> </tr> <tr> <td rowspan="2" style="text-align: center;">Depth</td> <td style="text-align: center;">16"</td> <td></td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">18"</td> <td></td> <td></td> <td></td> </tr> </table>		Select One!	Width					16"	18"	20"	Depth	16"				18"				\$198.00
	Select One!	Width																			
		16"	18"	20"																	
Depth	16"																				
	18"																				
<input type="checkbox"/> Nubuck Vinyl																					

\$

OPTIONAL ACCESSORIES - Additional Charges Apply

<input type="checkbox"/> Quad Cane Holder	\$86.40	<input type="checkbox"/> Rear Basket	\$72.00	<input type="checkbox"/> Stump Support*	\$450.00
<input type="checkbox"/> Cane Holder Tube	\$48.00	<input type="checkbox"/> Cup Holder	\$62.40	<input type="checkbox"/> Scooter Cover	\$117.60
<input type="checkbox"/> Crutch Holder Tube	\$84.00	<input type="checkbox"/> Pack N' Go	\$62.40	<input type="checkbox"/> Lap Belt	\$67.20
<input type="checkbox"/> Oxygen Tank Holder	\$132.00	<input type="checkbox"/> Pnuematic Rise Seat	\$98.00	<input type="checkbox"/> Flag	\$19.20
<input type="checkbox"/> Walker Holder	\$105.60				

* Not available on Standard Seat

\$

STANDARD FEATURES

18" W X 18" D Black Vinyl Stadium Seat	(1) 9" Flat Free Front Tire
(2) U1 Batteries	Side View Mirrors
90 Amp Controller	(2) 9" Flat Free Rear Tires
Battery Charger	Front Basket
Front and Rear Lighting Package	Height Adjustable Seat with Slide
Infinite Adjustable Delta Tiller	

SPECIFICATIONS

Model	Length	Width	Turning Radius	Weight Capacity Tested	Top Speed*	Range* (miles)
GC240 - 3 wheel	42.5"	23.5"	39"	220 lbs 300 lbs	5.0 mph 4.5 mph	19.5 miles 15.5 miles

* Top Speed and Range will vary due to rider weight, drive surface and drive slope.

WARRANTY

Metal Frame - Limited Lifetime	Battery Charger - 2 Year	Controller - 2 Year
Drive Train - 2 Years	Electrical Components** - 2 Year	** Does not include Batteries.

ASSIGNMENT OF A HCPCS CODE DOES NOT INDICATE COVERAGE OR PAYMENT FOR THE ITEM.

Total Cost

The intent of this form is as a Medicare Worksheet and not an Order Form. Prices listed are the Golden Technologies MSRP for each product and they are not the allowables for Medicare.