

Evaluated By: _____
 Mark For: _____
 Date: _____
 Account #: _____

Company: _____
 Address: _____
 Phone: _____

GC440 - COMPANION I SCOOTER - (4) Three Wheel K0807

GC440 COMPANION II - 4 Wheel Scooter
 350 lb Capacity with 18" X 16" Black Vinyl Stadium Seat **\$2,899.00**

STANDARD COLORS - Please Choose One!

Red **Blue**

OPTIONAL SEATING - Additional Charges Apply

| | | | | | | | | | | | | | | | | |
|---|--|---------------------------------|--------------|--|--|-----|-----|-----|--------------|-----|--|--|-----|--|--|--|
| Please Choose Color | Stadium Style Seat Sizes | Two-Tone Captain's Chair | | | | | | | | | | | | | | |
| <input type="checkbox"/> Charcoal Vinyl | <table border="1"> <tr> <td rowspan="2">Select One!</td> <td colspan="3">Width</td> </tr> <tr> <td>16"</td> <td>18"</td> <td>20"</td> </tr> <tr> <td rowspan="2">Depth</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>18"</td> <td></td> <td></td> </tr> </table> | Select One! | Width | | | 16" | 18" | 20" | Depth | 16" | | | 18" | | | <input type="checkbox"/> 18" W X 18" D High Back |
| Select One! | Width | | | | | | | | | | | | | | | |
| | 16" | 18" | 20" | | | | | | | | | | | | | |
| Depth | 16" | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | |
| <input type="checkbox"/> Nubuck Vinyl | | \$198.00 | | | | | | | | | | | | | | |

OPTIONAL ACCESSORIES - Additional Charges Apply

| | | |
|---|---|---|
| <input type="checkbox"/> Quad Cane Holder \$86.40 | <input type="checkbox"/> Rear Basket \$72.00 | <input type="checkbox"/> Stump Support* \$450.00 |
| <input type="checkbox"/> Cane Holder Tube \$48.00 | <input type="checkbox"/> Cup Holder \$62.40 | <input type="checkbox"/> Scooter Cover \$117.60 |
| <input type="checkbox"/> Crutch Holder Tube \$84.00 | <input type="checkbox"/> Pack N' Go \$62.40 | <input type="checkbox"/> Lap Belt \$67.20 |
| <input type="checkbox"/> Oxygen Tank Holder \$132.00 | <input type="checkbox"/> Pnuematic Rise Seat \$98.00 | <input type="checkbox"/> Flag \$19.20 |
| <input type="checkbox"/> Walker Holder \$105.60 | | |

** Not available on Standard Seat*

STANDARD FEATURES

| | |
|--|-----------------------------------|
| 18" W X 18" D Black Vinyl Stadium Seat | (2) 10" Flat Free Front Tires |
| (2) U1 Batteries | Side View Mirrors |
| 90 Amp Controller | (2) 10" Flat Free Rear Tires |
| Battery Charger | Front Basket |
| | Height Adjustable Seat with Slide |
| | Front and Rear Lighting Package |
| | Infinite Adjustable Delta Tiller |

SPECIFICATIONS

| Model | Length | Width | Turning Radius | Weight Capacity Tested | Top Speed | Range* (miles) |
|------------------------|--------|-------|----------------|------------------------|--------------------|--------------------------|
| GC440 - 4 wheel | 47.5" | 24" | 53.5" | 220 lbs 350 lbs | 5.0 mph 4.5 mph | 18.5 miles 15.0 miles |

* Battery range tested at 350 lb capacity. It will vary due to rider weight, drive surface and drive slope.

WARRANTY

| | | |
|---------------------------------------|---|---------------------------------------|
| Metal Frame - Limited Lifetime | Battery Charger - 2 Year | Controller - 2 Year |
| Drive Train - 2 Years | Electrical Components** - 2 Year | ** Does not include Batteries. |

ASSIGNMENT OF A HCPCS CODE DOES NOT INDICATE COVERAGE OR PAYMENT FOR THE ITEM.

The intent of this form is as a Medicare Worksheet and not an Order Form. Prices listed are the Golden Technologies MSRP for each product and they are not the allowables for Medicare.

| |
|-------------------|
| Totals |
| \$ |
| \$ |
| \$ |
| \$ |
| Total Cost |