

Evaluated By: _____
 Mark For: _____
 Date: _____
 Account #: _____

Company: _____
 Address: _____
 Phone: _____

GOLDEN COMPASS HD™ POWER BASE with Van / Pan Seat - Center-Wheel Drive K0824

GP620SS GOLDEN COMPASS HD POWER BASE with Van / Pan Seat - Center-Wheel Drive
 450 lb weight capacity **\$7,347.00**

STANDARD COLORS - Please Choose One!

Red Blue Black Green

BASIC EQUIPMENT PACKAGE FOR POWER WHEELCHAIRS

- Off - Board Battery Charger Adjustable Footplate Lap Belt (2) 6" Flat Free Front Casters
- 75 Amp Remote Controller Adjustable Armrests (2) 12" Flat Free Drive Wheels
- 20" W X 20" D Two-Tone Van / Pan Seat - Charcoal Vinyl (2) 6" Flat Free Rear Casters

***Note:** All Golden Power Wheelchairs are shipped with Programable Controllers.
 Non-Programmable Controllers are available by request only!

SPECIFICATIONS

Model	Length	Width	Turning Radius	Top Speed	Range* (miles)
GP620SS - Center Wheel	39.8"	25"	20"	4 mph	12

* Battery range will vary due to rider weight, drive surface and drive slope

WARRANTY

- Metal Frame - 5 Year Battery Charger - 1 Year Controller - 1 Year
- Drive Train - 1 Year Electrical Components** - 1 Year ** Does not include Batteries.
- Plastic Parts (excludes shroud), Rubber Parts (excludes tires), and Bearings - 1 Year

SEATING OPTIONS - Additional Charges May Apply

- Two-Tone Van Pan Seat Proflex Seating Frame **\$500.00**
 - Gray/Charcoal Vinyl
 - 22" W X 20" D **\$198.00**
 - 24" W X 20" D **\$198.00**
- | Select One! | Width | | |
|-------------|-------|-----|-----|
| | 18" | 20" | 22" |
| Depth | 18" | | |
| | 20" | | |

OPTIONS and ACCESSORIES - Additional Charges May Apply

- (2) 12 Volt/22NF Batteries E2361 **\$500.00** Elevating Leg Rests E0990 **\$391.20**
- Oxygen Tank Holder E2208 **\$132.00** Walker Holder K0108 **\$110.40**
- Quad Cane Holder E2207 **\$116.40** Stump Support E1020 **\$333.60**
- Cane Holder Tube E2207 **\$48.00** Powerchair Cover **\$117.60**
- Crutch Holder Tube E2207 **\$84.00** Pack N' Go **\$62.40** Flag **\$19.20**
- Swing Away Footrests K0052* **\$225.00** * Replacement Only Rear Basket **\$72.00**

ASSIGNMENT OF A HCPCS CODE DOES NOT INDICATE COVERAGE OR PAYMENT FOR THE ITEM.

The intent of this form is as a Medicare Worksheet and not an Order Form. Prices listed are the Golden Technologies MSRP for each product and they are not the allowables for Medicare.

Totals
\$
\$
\$
\$
Total Cost