



LM502 LIFT FRAME PARTS LIST

Part Number	Description	Qty. Per Unit
BB5506	1/4" LOCK WASHER	17
BB5514	HAIR PIN	12
BB5530	5/16 X 1 LEG LEVELER	13
BB5535	1/8 NYLON PLASTIC WASHER	6
BB5548	5/16" X 3" T AP BOLT GRADE 5 (to attach lift frame to chair)	4
BB5549	3/8 NYLON PLASTIC WASHER	4
BB5554	5/16" NYLON LOCK NUT	6
BB5583	1 1/4" X 1 1/4" END CAP	12
BB5584	1/2" X 1/2" SHOULDER BOLT	2
BB5585	1/2" X 3 1/4" CLEVIS PIN	2
BB5586	1/2" X 2 1/2" SHOULDER BOLT	2
BB5587	3/8" X 1 1/4" SHOULDER BOLT	1
BB5588	3/8"-16 NYLON LOCK NUT	4
BB5589	3/8" X 1 3/4" SHOULDER BOLT	3
BB5590	3/8" X 2 1/2" SHOULDER BOLT	2
BB5591	5/16" X 1 1/2" HEX CAP SCREW (to attach back to mechanis	2
BB5592	5/16" X 2" HEX CAP SCREW 9 to (to attach back to mechanis	4
BB5593	NYLON WASHER .750 X .502 X .1875	8
BB5594	NYLON WASHER 1.25 X .502 X .0625	2
BB5596	1/2" X 4" CLEVIS PIN	4
BB5597	1/2" X 2 1/2" CLEVIS PIN	2
BB5598	1/2" X 2 1/4" CLEVIS PIN	4
BB5599	5/8" X 1/2" X 1/2" BUSHING	28
BB5599KIT	BUSHING KIT FOR 502	1
GMOKI-28	26-28 WIDE MOTOR	2
GMOKICB	26-28 WIDE CONTROL BOX	1
GMOKI1502C	SINGLE MOTOR 502	1
GM502Y	Y CONNECTOR	1
LIFT-502	COMPLETE 502 LIFT FRAME ASSEMBLY (Recline bar not included)	1

LM502	502 LIFT FRAME	1
LM5021	TOP PIECE 502 LIFT FRAME	1
LM5022	STRAIGHT BAR 502 LIFT FRAME	1
LM5023	H-BAR 502 LIFT FRAME	1
LM5024	BASE PIECE 502 LIFT FRAME	1
RS502	RECLINE STRAP 502 LIFT FRAME	4
ZK756HC	HAND CONTROL FOR 756/757/502	1

RECLINE BAR FOR 502 LIFT FRAME

RB502	RECLINE BAR FOR 502 CHAIR	1
-------	---------------------------	---

RECLINE MECHANISM FOR 502 LIFT FRAME

RM502BR	BACK RECLINE MECHANISM 502	1
RM502LH	RECLINE MECHANISM LEFT SIDE 502	1
RM502RH	RECLINE MECHANISM RIGHT SIDE 502	1
	LEFT AND RIGHT IS DETERMINED BY SITTING IN CHAIR	

TRANSFORMER FOR 502 SERIES LIFT FRAME

GM1501	TRANSFORMER-D	2
--------	---------------	---

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING



LM5100 LIFT FRAME PARTS LIST

Part Number	Description	Qty. Per Unit
BB5003	1/16" NYLON PLASTIC WASHER	2
BB5512	3/8" x 1-5/8" CLEVIS PIN	1
BB5514	HAIR PIN	12
BB5515	1 INCH SQUARE END CAP	8
BB5518	8 INCH TYWRAP	2
BB5530	5/16" x 1" LEG LEVELER	4
BB5531	3/8" x 2-1/2" CLEVIS PIN	2
BB5535	1/8" NYLON PLASTIC WASHER	8
BB5549	3/8" NYLON PLASTIC WASHER	2
BB5566	3/8" x 1-3/4" CLEVIS PIN	4
BB5567	3/8" x 2-5/8" CLEVIS PIN	1
BB5579	3/8" x 2-7/8" CLEVIS PIN	4
GM1500	DC INLINE MOTOR	1
GM1500	Assy DC MOTOR COMPLETE	1
LIFT -5100	COMPLETE 5100 LIFT FRA (Recline bar not included)	1
LM5100	5100 LIFT FRAME	1
LM5101	TOP PIECE 5100 LIFT FRAME	1
LM5102	U BAR 5100 SERIES LIFT FRAME	1
LM5103	STRAIGHT BAR 5100 SERIES LIFT FRAME	2
LM5104	BASE PIECE 5100 LIFT FRAME	1
LM5100BP	BRIDGE PLATE 5100 SERIES LIFT	2
RS5100	RECLINE STRAP 5100 / 6111 SERIES LIFT FRAME	2
ZKOKI-DQ	QUICK RELEASE DC PENDANT	1

RECLINE BARS FOR 5100 LIFT FRAME

RB753-1	RECLINE BAR PETIT / 753	1
RB4000-C	RECLINE BAR CAPRI	1

RB5100-M

RECLINE BAR MEDIUM / SMALL

1

RM621L

3 WAY MECHANISM LEFT SIDE

1

RM621R

3 WAY MECHANISM RIGHT SIDE

1

RM8300L

2 WAY MECHANISM LEFT SIDE

1

RM8300R

2 WAY MECHANISM RIGHT SIDE

1

LEFT AND RIGHT IS DETERMINED BY SITTING IN CHAIR

GM1501

TRANSFORMER FOR 5100 LIFT FRAME

TRANSFORMER-D

1

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

LM5120 LIFT FRAME PARTS LIST

Part Number	Description	Qty. Per Unit
BB5003	1/16" NYLON PLASTIC WASHER	2
BB5512	3/8" x 1-5/8" CLEVIS PIN	1
BB5514	HAIR PIN	12
BB5515	1 INCH SQUARE END CAP	8
BB5518	8 INCH TYWRAP	2
BB5530	5/16" x 1" LEG LEVELER	4
BB5531	3/8" x 2-1/2" CLEVIS PIN	2
BB5535	1/8" NYLON PLASTIC WASHER	8
BB5549	3/8" NYLON PLASTIC WASHER	2
BB5566	3/8" x 1-3/4" CLEVIS PIN	4
BB5567	3/8" x 2-5/8" CLEVIS PIN	1
BB5579	3/8" x 2-7/8" CLEVIS PIN	4
GM1500	DC INLINE MOTOR	1
GM1500	Assy DC MOTOR COMPLETE	1
LIFT -5120	COMPLETE 5120 LIFT FRA (Recline bar not included)	1
LM5120	5120 LIFT FRAME	1
LM5121	TOP PIECE 5120 LIFT FRAME	1
LM5102	U BAR 5100 SERIES LIFT FRAME	1
LM5103	STRAIGHT BAR 5100 SERIES LIFT FRAME	2
LM5124	BASE PIECE 5100 LIFT FRAME	1
LM5100BP	BRIDGE PLATE 5100 SERIES LIFT	2
RS5100	RECLINE STRAP 5100 / 6111 SERIES LIFT FRAME	2
ZKOKI-DQ	QUICK RELEASE DC PENDANT	1
<h3>RECLINE BARS FOR 5120 LIFT FRAME</h3>		
RB5100-M	RECLINE BAR MEDIUM / SMALL	1
RB5100-L	RECLINE BAR LARGE	1

RECLINE MECHANISMS FOR 5120 LIFT FRAME

RM621L
RM621R

3 WAY MECHANISM LEFT SIDE 1
3 WAY MECHANISM RIGHT SIDE 1
LEFT AND RIGHT IS DETERMINED BY SITTING IN CHAIR

TRANSFORMER FOR 5120 LIFT FRAME

GM1501

TRANSFORMER-D 1

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING



LM6111 LIFT FRAME PARTS LIST

Part Number	Description	Qty. Per Unit
BB5003	1/16" NYLON PLASTIC WASHER	4
BB5512	3/8" x 1-5/8" CLEVIS PIN	2
BB5514	HAIR PIN	15
BB5515	1 INCH SQUARE END CAP	8
BB5518	8 INCH TYWRAP	2
BB5530	5/16" X 1" LEG LEVELER	6
BB5531	3/8" x 2-1/2" CLEVIS PIN	4
BB5535	1/8" NYLON PLASTIC WASHER	8
BB5549	3/8" NYLON PLASTIC WASHER	4
BB5566	3/8" x 1-3/4" CLEVIS PIN	4
BB5567	3/8" x 2-5/8" CLEVIS PIN	2
BB5582	3/8" x 3" CLEVIS PIN	4
GMOKI-28W	26" / 28" WIDE DC MOTOR	2
GMOKI-CB	26"-28" WIDE CONTROL	1
GM1502-EXT	(Extension wire low voltage)	1
LIFT -6111	COMPLETE 6111 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM6111	6111 LIFT FRAME	1
LM61111	TOP PIECE 6111 LIFT FRAME	1
LM61112	STRAIGHT BAR 6111 SERIES LIFT FRAME	2
LM61113	H BAR 6111 SERIES LIFT FRAME	1
LM61114	BASE PIECE 6111 LIFT FRAME	1
RS5100	RECLINE STRAP 5100 / 6111 SERIES LIFT FRAME	4
ZKOKI-DQ	QUICK RELEASE DC PENDANT	1

RECLINE BARS FOR 6111 LIFT FRAME

RB6111	DUAL MOTOR RECLINE BAR	1
RB6111-28	DUAL MOTOR RECLINE BAR 28 WIDE	1

RECLINE MECHANISMS FOR 6111 LIFT FRAME

RM621HD-L 3 WAY MECHANISM LEFT SIDE HEAVY DUTY	1
RM621HD-R 3 WAY MECHANISM RIGHT SIDE HEAVY DUTY	1
LEFT AND RIGHT IS DETERMINED BY SITTING IN CHAIR	

TRANSFORMER FOR 6111 LIFT FRAME

GM1501 TRANSFORMER-D	1
----------------------	---

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

LM756 LIFT FRAME PARTS LIST

Part Number	Description	Qty. Per Unit
BB5003	1/16" NYLON PLASTIC WASHER	4
BB5512	3/8" x 1-5/8" CLEVIS PIN	3
BB5514	HAIR PIN	15
BB5515	1 INCH SQUARE END CAP	8
BB5518	8 INCH TYWRAP	2
BB5530	5/16" x 1" LEG LEVELER	4
BB5531	3/8" x 2-1/2" CLEVIS PIN	1
BB5535	1/8" NYLON PLASTIC WASHER	8
BB5940	13/32" NYLON PLASTIC WASHER	2
BB5549	3/8" NYLON PLASTIC WASHER	8
BB5554	5/16" LOCK NUT	2
BB5566	3/8" x 1-3/4" CLEVIS PIN	4
BB5572	1-1/4" NYLON WASHER	2
BB5576	3/8" x 5/8" SHOUDER	2
BB5577	3/8" x 2-3/4" CLEVIS PIN	2
BB5579	3/8" x 2-7/8" CLEVIS PIN	4
GM756	756 LIFT MOTOR	2
LIFT -756	COMPLETE 756 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM756	756 LIFT FRAME	1
LM7561	TOP PIECE 756 LIFT FRAME	1
LM7562	STRAIGHT BAR 756 SERIES LIFT FRAME	2
LM7563	H BAR 756 SERIES LIFT FRAME	1
LM7564	BASE PIECE 756 LIFT FRAME	1
LM7565	BRIDGE PLATE 756 SERIES LIFT	2
RS756	RECLINE STRAP 756 SERIES LIFT FRAME	2
ZK756HC	756 HAND CONTROL	1

RECLINE BAR FOR 756 LIFT FRAME

RB756	RECLINE BAR 756 SERIES	1
-------	------------------------	---

RECLINE MECHANISMS FOR 756 LIFT FRAME

RM756BM	RECLINE MECHANISM BACK MOUNT 756 SERIES	2
RM756BR	BACKRECLINE MECHANISM 756 SERIES	1
RM756L	RECLINE MECHANISM LEFT SIDE 756 SERIES	1
RM756R	RECLINE MECHANISM RIGHT SIDE 756 SERIES	1

LEFT AND RIGHT IS DETERMINED BY SITTING IN CHAIR

TRANSFORMER FOR 756 SERIES LIFT FRAME

GM756CB	CONTROL BOX 756 SERIES	1
---------	------------------------	---

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

756 MC PARTS LIST
ELECTRICAL SYSTEM FOR 756MC

Part Number	Description	Qty. Per Unit
GM756L	OKIN LIFT MOTOR MAXI COMFORT	1
GM756R	OKIN BACK RECLINE MOTOR MAXI COMFORT	1
ZK756HCL	HAND CONTROL FOR 756/502/505	1
GM1501	TRANSFORMER-D	1

EXT100-DQ	OKIN HAND CONTROL EXTENSION	1
GM756-ASSY	DC MOTORS COMPLETE FOR MAXI COMFORT	1
BB5003	1/16 NYLON PLASTIC WASHER	2
BB5549	5/16 NYLON PLASTIC WASHER	2
BB5531	3/8 X 2 1/2 CLEVIS PIN	1
BB5512	3/8 X 1 5/8 CLEVIS PIN	3
BB5514	HAIR PIN	4

See LM5120 lift motor assembly and LM5120 lift frame parts list for assembly and parts identification

LIFT FRAME FOR 756MC

LIFT-5120	COMPLETE 5120 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM5120	5120 LIFT FRAME	1

See LM5120 lift frame and LM5120 lift frame parts list for assembly and parts identification

RECLINE BAR FOR 756MC

RB5100-M	RECLINE BAR MEDIUM	1
----------	--------------------	---

See LM5120 lift frame for assembly and parts identification

RECLINE MECHANISM FOR 756MC

RM620L	INFINITE-3 WAY MECH.-LEFT	1
RM620R	INFINITE-3 WAY MECH.-RIGHT	1

RM756FB	RM756 FRONT BRACKET	1
RM756RMB	RM756 REAR MOUNTING BRACKET -YOKE	1
RM756SB	RM756 SPACER BAR	1
BB5608	620 MECH STOP	2
BB5603	1/4-20 NYLON LOCK NUT FLANGE	8
BB5602	1/4X20X1.5 HEX HEAD FLANGE BOLT	4
BB5605	1/4 X 20 X 1 FLAT HEAD BOLT	4
BB5607	1/4 X 20 X 3/4 HEX CAP SCREW	4
BBWASH	5/16 LOCK WASHER	2
MS1.5WL	1/4 X 1 1/2 WHIZ LOCK BOLT	4
MS2.5WL	1/4 X 2 1/2 WHIZ LOCK BOLT	2
MS3.0TP	TAP BOLT (GRADE 5-3 INCH See Maxi Comfort recline mechanism for assembly	2

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

505 SMALL PARTS LIST
ELECTRICAL SYSTEM FOR 505 SMALL

Part Number	Description	Qty. Per Unit
GM756L	OKIN LIFT MOTOR MAXI COMFORT	1
GM756R	OKIN BACK RECLINE MOTOR MAXI COMFORT	1
ZK756HCL	HAND CONTROL FOR 756/502/505	1
GM1501	TRANSFORMER-D	1
EXT100-DQ	OKIN HAND CONTROL EXTENSION	1
GM756-ASSY	DC MOTORS COMPLETE FOR MAXI COMFORT	1
BB5003	1/16 NYLON PLASTIC WASHER	2
BB5549	5/16 NYLON PLASTIC WASHER	2
BB5531	3/8 X 2 1/2 CLEVIS PIN	1
BB5512	3/8 X 1 5/8 CLEVIS PIN	3
BB5514	HAIR PIN	4

See LM5100 lift motor assembly and LM5100 lift frame parts list for assembly and parts identification

LIFT FRAME FOR 505 SMALL

LIFT-5100	COMPLETE 5100 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM5100	5100 LIFT FRAME	1

See LM5100 lift frame and LM5100 lift frame parts list for assembly and parts identification

RECLINE BAR FOR 505 SMALL

RB5100-S	RECLINE BAR SMALL	1
----------	-------------------	---

See LM5100 lift frame for assembly and parts identification

RECLINE MECHANISM FOR 505 SMALL

RM620L	INFINITE-3 WAY MECH.-LEFT	1
RM620R	INFINITE-3 WAY MECH.-RIGHT	1

RM758SFB	RM758S FRONT BRACKET	1
RM758SRMB	RM758S REAR MOUNTING BRACKET-YOKE	1

RM758SSB	RM758S SPACER BAR	1
BB5608	620 MECH STOP	2
BB5603	1/4-20 NYLON LOCK NUT FLANGE	8
BB5602	1/4X20X1.5 HEX HEAD FLANGE BOLT	4
BB5605	1/4 X 20 X 1 FLAT HEAD BOLT	4
BB5607	1/4 X 20 X 3/4 HEX CAP SCREW	4
BBWASH	5/16 LOCK WASHER	2
MS1.5WL	1/4 X 1 1/2 WHIZ LOCK BOLT	4
MS2.5WL	1/4 X 2 1/2 WHIZ LOCK BOLT	2

MS3.0TP	TAP BOLT (GRADE 5-3 INCH	2
	See Maxi Comfort recline mechanism for assembly	

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

505 MEDIUM PARTS LIST
ELECTRICAL SYSTEM FOR 505 MEDIUM

Part Number	Description	Qty. Per Unit
GM756L	OKIN LIFT MOTOR MAXI COMFORT	1
GM756R	OKIN BACK RECLINE MOTOR MAXI COMFORT	1
ZK756HCL	HAND CONTROL FOR 756/502/505	1
GM1501	TRANSFORMER-D	1

EXT100-DQ	OKIN HAND CONTROL EXTENSION	1
GM756-ASSY	DC MOTORS COMPLETE FOR MAXI COMFORT	1
BB5003	1/16 NYLON PLASTIC WASHER	2
BB5549	5/16 NYLON PLASTIC WASHER	2
BB5531	3/8 X 2 1/2 CLEVIS PIN	1
BB5512	3/8 X 1 5/8 CLEVIS PIN	3
BB5514	HAIR PIN	4

See LM5100 lift motor assembly and LM5100 lift frame parts list for assembly and parts identification

LIFT FRAME FOR 505 MEDIUM

LIFT-5100	COMPLETE 5100 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM5100	5100 LIFT FRAME	1

See LM5100 lift frame and LM5100 lift frame parts list for assembly and parts identification

RECLINE BAR FOR 505 MEDIUM

RB5100-M	RECLINE BAR MEDIUM	1
----------	--------------------	---

See LM5100 lift frame for assembly and parts identification

RECLINE MECHANISM FOR 505 MEDIUM

RM620L	INFINITE-3 WAY MECH.-LEFT	1
RM620R	INFINITE-3 WAY MECH.-RIGHT	1

RM756FB	RM756 FRONT BRACKET	1
RM756RMB	RM756 REAR MOUNTING BRACKET -YOKE	1

RM756SB	RM756 SPACER BAR	1
BB5608	620 MECH STOP	2
BB5603	1/4-20 NYLON LOCK NUT FLANGE	8
BB5602	1/4X20X1.5 HEX HEAD FLANGE BOLT	4
BB5605	1/4 X 20 X 1 FLAT HEAD BOLT	4
BB5607	1/4 X 20 X 3/4 HEX CAP SCREW	4
BBWASH	5/16 LOCK WASHER	2
MS1.5WL	1/4 X 1 1/2 WHIZ LOCK BOLT	4
MS2.5WL	1/4 X 2 1/2 WHIZ LOCK BOLT	2

MS3.0TP	TAP BOLT (GRADE 5-3 INCH	2
	See Maxi Comfort recline mechanism for assembly	

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING

505 LARGE PARTS LIST
ELECTRICAL SYSTEM FOR 505 LARGE

Part Number	Description	Qty. Per Unit
GM756L	OKIN LIFT MOTOR MAXI COMFORT	1
GM756R	OKIN BACK RECLINE MOTOR MAXI COMFORT	1
ZK756HCL	HAND CONTROL FOR 756/502/505	1
GM1501	TRANSFORMER-D	1
EXT100-DQ	OKIN HAND CONTROL EXTENSION	1
GM756-ASSY	DC MOTORS COMPLETE FOR MAXI COMFORT	1
BB5003	1/16 NYLON PLASTIC WASHER	2
BB5549	5/16 NYLON PLASTIC WASHER	2
BB5531	3/8 X 2 1/2 CLEVIS PIN	1
BB5512	3/8 X 1 5/8 CLEVIS PIN	3
BB5514	HAIR PIN	4

See LM5120 lift motor assembly and LM5120 lift frame parts list for assembly and parts identification

LIFT FRAME FOR 505 LARGE

LIFT-5120	COMPLETE 5120 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM5120	5120 LIFT FRAME	1

See LM5120 lift frame and LM5120 lift frame parts list for assembly and parts identification

RECLINE BAR FOR 505 LARGE

RB5100-L	RECLINE BAR LARGE	1
----------	-------------------	---

See LM5120 lift frame for assembly and parts identification

RECLINE MECHANISM FOR 505 LARGE

RM620L	INFINITE-3 WAY MECH.-LEFT	1
RM620R	INFINITE-3 WAY MECH.-RIGHT	1

RM758LFB	RM758L FRONT BRACKET	1
RM758LRMB	RM758L REAR MOUNTING BRACKET -YOKE	1

RM758LSB	RM758L SPACER BAR	1
BB5608	620 MECH STOP	2
BB5603	1/4-20 NYLON LOCK NUT FLANGE	8
BB5602	1/4X20X1.5 HEX HEAD FLANGE BOLT	4
BB5605	1/4 X 20 X 1 FLAT HEAD BOLT	4
BB5607	1/4 X 20 X 3/4 HEX CAP SCREW	4
BBWASH	5/16 LOCK WASHER	2
MS1.5WL	1/4 X 1 1/2 WHIZ LOCK BOLT	4
MS2.5WL	1/4 X 2 1/2 WHIZ LOCK BOLT	2

MS3.0TP	TAP BOLT (GRADE 5-3 INCH	2
	See Maxi Comfort recline mechanism for assembly	

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD P ART NUMBERS AND PRICING

505 TALL PARTS LIST
ELECTRICAL SYSTEM FOR 505 TALL

Part Number	Description	Qty. Per Unit
GM756L	OKIN LIFT MOTOR MAXI COMFORT	1
GM756R	OKIN BACK RECLINE MOTOR MAXI COMFORT	1
ZK756HCL	HAND CONTROL FOR 756/502/505	1
GM1501	TRANSFORMER-D	1
EXT100-DQ	OKIN HAND CONTROL EXTENSION	1
GM756-ASSY	DC MOTORS COMPLETE FOR MAXI COMFORT	1
BB5003	1/16 NYLON PLASTIC WASHER	2
BB5549	5/16 NYLON PLASTIC WASHER	2
BB5531	3/8 X 2 1/2 CLEVIS PIN	1
BB5512	3/8 X 1 5/8 CLEVIS PIN	3
BB5514	HAIR PIN	4

See LM5120 lift motor assembly and LM5120 lift frame parts list for assembly and parts identification

LIFT FRAME FOR 505 TALL

LIFT-5120	COMPLETE 5120 LIFT FRAME ASSEMBLY (Recline bar not included)	1
LM5120	5120 LIFT FRAME	1

See LM5120 lift frame and LM5120 lift frame parts list for assembly and parts identification

RECLINE BAR FOR 505 TALL

RB5100-L	RECLINE BAR LARGE	1
----------	-------------------	---

See LM5120 lift frame for assembly and parts identification

RECLINE MECHANISM FOR 505 TALL

RM620L	INFINITE-3 WAY MECH.-LEFT	1
RM620R	INFINITE-3 WAY MECH.-RIGHT	1

RM758LFB	RM758L FRONT BRACKET	1
RM758LRMB	RM758L REAR MOUNTING BRACKET-YOKE	1

RM758LSB	RM758L SPACER BAR	1
BB5608	620 MECH STOP	2
BB5603	1/4-20 NYLON LOCK NUT FLANGE	8
BB5602	1/4X20X1.5 HEX HEAD FLANGE BOLT	4
BB5605	1/4 X 20 X 1 FLAT HEAD BOLT	4
BB5607	1/4 X 20 X 3/4 HEX CAP SCREW	4
BBWASH	5/16 LOCK WASHER	2
MS1.5WL	1/4 X 1 1/2 WHIZ LOCK BOLT	4
MS2.5WL	1/4 X 2 1/2 WHIZ LOCK BOLT	2

MS3.0TP	TAP BOLT (GRADE 5-3 INCH	2
	See Maxi Comfort recline mechanism for assembly	

CALL TECHNICAL SERVICE FOR FABRIC, FOAM, SPRING AND WOOD PART NUMBERS AND PRICING