

P O W E R C H A I R



- Seat Colors
- Charcoal/Gray Vinyl
 - Coffee/Sand Vinyl
- Body Colors
- Candy Apple Red
 - Blue

Roll

New!

golden COMPASS[™] SPORT

The Golden Compass family now includes the all-new Compass Sport! Featuring the same center-wheel drive performance as the original Compass, the Compass Sport features a sleek new design to the shroud and a single post seat mount. The Compass Sport comes standard in Red and Blue.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

Golden Compass™ Sport Specifications

Model Number	GP605CC	GP605SS
Medicare Code	K0823	K0822
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Center	Center
Maximum Speed	4 mph	4 mph
Operating Range ³	15 Miles	15 Miles
Ground Clearance to transaxle/motor	3.25"	3.25"
Turning Radius	19.5"	19.5"
Type Batteries	2-U1	2-U1
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	Dynamic Shark 60 Amp Remote	Dynamic Shark 60 Amp Remote
Dimensions:		
Length	39"	39"
Width	24"	24"
Height (ground to top of back of seat)	40" - 42"	40" - 42"
Height (ground to top of back of headrest)	47" - 52.5"	47" - 52.5"
Ground to Top of Seat	21" - 23"	19" - 21"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Weight of Unit: (assembled)		
Base	92 lbs	92 lbs
Seat w/Arms	43 lbs	40 lbs
Batteries	22 lbs ea	22 lbs ea
Tires:		
Rear	10"	10"
Casters	6"	6"
Standard Seat:		
Slide Seat	Captain's Highback	Captain's Van/Pan
Back Height - no headrest	N/A	N/A
Back Height with headrest	21"	23"
Width x Depth	16" X 16", 18" X 18", 20" X 20"	18" X 18", 20" X 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
Optional Seats:		
Back Height - no headrest	Captain's Highback	Captain's Van/Pan
Back Height with headrest	21"	23"
Width	16", 18", 20" 22"	16", 18", 20" 22"
Depth	16", 18" or 20"	16", 18" or 20"
Color	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray
Optional Seats 2:		
Back Height - no headrest	N/A	Pro-Flex Seat
Back Height with headrest	N/A	N/A
Width	N/A	16", 18", 20", or 22"
Depth	N/A	16", 18" or 20"
Color	N/A	Black
Seat Options:		
Pneumatic Rise Seat Option	N/A	N/A
Electric Rise Seat Option	N/A	N/A
Footrest	Yes	Yes
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 - 35 degrees	0 - 35 degrees
Warranty:		
Frame	Lifetime	Lifetime
Drive Train	13 Months	13 Months
Electronics ²	13 Months	13 Months



Compass Sport GP605SS shown

2 exciting colors



Candy Apple Red



Blue

Available Accessories:

- Crutch Holder Tube
- Stump Support
- Cane Holder Tube
- Elevating Leg Rests
- Oxygen Tank Holder
- Power Chair Cover
- Quad Cane Holder
- Rear Basket
- Walker Holder
- Pack N' Go
- Swing-away Desk Arm
- Van/pan Seat
- Swing-away Footrests

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

³ Range listed tested at 200 lbs.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.