

LIGHTWEIGHT PORTABLE SCOOTER



Body Colors



Red



Blue



Sahara

**New tire & wheel design!**

# Ride

**BUZZ**around<sup>®</sup>  
*Lite*

3 wheel  
GB-106

The popular Buzzaround Lite scooter has been improved! This compact scooter now comes standard with a Delta tiller with wraparound handle and a larger, more comfortable seat. The Buzzaround Lite still features our easy "Pop! Pull! Push!" wireless disassembly for transport. Change the look of your scooter with our three sets of snap-on shroud panels in Red, Blue and Sahara, all while enjoying its tight indoor maneuverability and smooth outdoor use.



Rise.

Roll.

Rest.

Ride.

**GOLDEN**  
technologies

YOUR LIFE IN MOTION

### 3 wheel Specifications

<b>Model</b>	<b>GB106</b>
Medicare Code	K0800
Weight Capacity	250 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range <sup>1</sup>	8 miles
Ground Clearance to center deck	1.75"
Ground Clearance to transaxle/motor	1.5"
Turning Radius	31"
Type Batteries	2-12AH
Colors	Red, Blue, Sahara
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	45 Amp
<b>Dimensions:</b> <sup>3</sup>	
Length	36.5"
Width	21.25"
Height (ground to top of back of seat)	30.375" - 32.375"
Ground to Top of Seat	19.25" - 21.25"
Ground to Top of Deck	4.25"
Deck to Top of Seat	15" - 17"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	28.625"
<b>Weight of Unit: (assembled)</b>	96 lbs
Front Half	29 lbs
Rear Half	26 lbs
Seat w/Arms	21 lbs
Batteries	19 lbs/pack
<b>Tires:</b>	Solid
Front	8"
Rear	8"
Rear Anti-Tippers	2"
<b>Standard Seat:</b>	Stadium Style
Back Height - no headrest	16"
Width x Depth	17" X 16"
Color	Black Vinyl
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years

<sup>1</sup>. Battery range will vary due to rider weight, drive surface and terrain.

<sup>2</sup>. Electronics warranty excludes batteries.

<sup>3</sup>. There is a +1"/-1" tolerance on all dimensions.



Exclusive interlocking teeth provide the tightest, most dependable battery connection for lasting power.

### Wireless, quick disconnect in three easy steps:



**Step 1:** Pop off the seat!



**Step 2:** Pull off the battery pack!




**Step 3:** Push apart the frame!



Quickly disassembles into five, lightweight pieces for easy transport!

### Features

- Exclusive interlocking teeth for fully connectorless assembly/disassembly
- 250 lb. weight capacity
- Interchangeable color panels included with each unit—FREE!
- Auto grade cut pile carpet
- Tiller mounted basket



New control panel design with Delta tiller standard