


MEDTRADE EDITION

M A X I C O M F O R T <sup>TM</sup> S E R I E S



\* upcharge applies

**Pub  
PR713**

The Pub chair is a traditional style two pillow waterfall back recliner featuring decorative brass nail head trim on both the arms and wings. The supple waterfall, zipper back design allows for customizing comfort. This dual motor chair is part of the MaxiComfort series. Recline in style!



# Pub Chair

	Pub
Model Number	PR713
Size	Medium
Weight Capacity	375 Lb.
Recline Positions	MaxiComfort
Standard Fabrics	3
Fabric Options	ALL
Chaise Seat	YES
Coil Springs	YES
Back Type	Pillow
Head Rest Cover	NO
Arm Covers	NO
Battery Back Up	YES
Overall Width	35"
Overall Height	45"
Width Between Arms	21"
Floor to Top of Seat	20"
Seat Depth	18.5"
Seat to Top of Back	26.5"
Seat to Top of Arm	4.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	26"
Back Angle Closed	100°
Back Angle Reclined	170°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	26"
Extension reclined from Back of Seat to end of Foot Rest	43.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	NO
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	NO
Extra Pocket Option Available	NO
Heat & Massage Option Available	NO
Gel Memory Foam Option Available	NO
Power Pillow Option Available	NO
AutoFlex Option Available	NO
Foot Rest Extension Option Available	NO
Straight Lift Option Available	NO
User Height (suggested)	5'4" - 5'10"



Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



## Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



## Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with  
UL 73 and  
CSA C22.2 No. 68

