

Buzzaround EX

Bumper Conversion Kit & Battery Docking Station



BUZZ[™]
AROUND
EXTREME

 **GOLDEN**



Bumper Conversion Kit & Battery Docking Station

Features

- **Control Panel:**
Easy to see LED battery gauge; easy accessible speed control; key switch and horn button.
- **Delta Tiller handle:**
Offers comfort and easy right or left hand use.
- **Headlight On/Off Switch:**
Located on the right side of the tiller for easy access.
- **Front Basket:**
Easily store your items.
- **LED Headlight:**
Angle-adjustable headlight, easier to see and be seen.
- **Tiller Adjustment and Lock Knobs:**
Lower and lock the tiller in place for more compact storage without removing the battery packs.
- **Optional Battery Docking Station**



Full Suspension:
Front and rear springs offer an all-day comfort ride.

BUMPER CONVERSION KIT FOR MODEL GB118 ONLY.



Front Bumper and Side Rails:
Protect the scooter completely from impacts.



Large Rear Bumper:
Offers added protection from mishaps. The caster wheels offer ease of storage by allowing the scooter to be rolled into place while standing vertically.

Storage:
Offers vertical storage with a compact footprint of 22 inches wide by 23 inches deep.

New!

Battery Docking Station for Models GB118 & GB148

The new Battery Docking Station is portable and charges the interlocking U1 batteries. Perfect for keeping a set of charged batteries on hand.

