

LIGHTWEIGHT PORTABLE SCOOTER



Comfort Spring System!



Seat Color
Black Vinyl

Body Colors
Blue
Red



3 wheel GB117S
4 wheel GB147S

The next generation of the popular Buzzaround scooter series is now available with our comfort spring system! The Buzzaround XLS offers more standard features to make your riding experience better than ever. Features a wraparound Delta tiller, a new LED battery gauge and a new LED angle-adjustable headlight. The new Buzzaround XLS is also longer to offer more foot room, with a great turning radius. The Buzzaround XLS disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!





LIGHTWEIGHT PORTABLE SCOOTER

Item	GB117S	GB147S
Description	XLS	XLS
Medicare Code	N/A	N/A
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to ¹	8 miles	8 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-12AH	2-12AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp	50 Amp
Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23"
Ground to Top of Deck	4.75"	4.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"
Weight of Unit: (assembled)³		
Front Half	105 lbs.	113 lbs.
Rear Half	35 lbs.	43 lbs.
Seat w/Arms	28 lbs.	28 lbs.
Batteries	23 lbs.	23 lbs.
Tires:		
Front	19 lbs. pack	19 lbs. pack
Rear	Solid	Solid
Rear Anti-Tip Wheels	8"	8"
Standard Seat		
Seat Back Height - no headrest	2"	2"
Width x Depth	Stadium Style	Stadium Style
Color	16"	16"
Optional Seat:		
Back Height - no headrest	18" x 16"	18" x 16"
Width x Depth	Black Vinyl	Black Vinyl
Color	16"	16"
Warranty:		
Frame	20" x 18"	20" x 18"
Drive Train	Black Vinyl	Black Vinyl
Electronics ²	Lifetime Warranty	Lifetime Warranty
Batteries	2 Year Warranty	2 Year Warranty
	2 Year Warranty	2 Year Warranty
	1 Year Warranty	1 Year Warranty

¹. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.
². Electronics warranty excludes batteries.



Comfort Spring System


Wireless disassembly in 4 easy steps!



Step 1: Pop off the seat!



Step 2: Pull out the battery pack!



Step 4: Pull forward on the drivetrain release lever to separate the front and rear frame sections!

Step 3: Fold down tiller and lock in place!



Disassembles into 5 pieces!

4 wheel model shown

Features

- 300 lb. weight capacity
- Angle-adjustable LED headlight
- Tiller mounted basket
- Comfort Spring System

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

