

S I G N A T U R E S E R I E S



**ALL NEW
USB Port
Standard!**



Anchor



Bittersweet



Sandstorm



Coffee Bean*



Buckskin*

* upcharge applies



PR405-MLA Orion
with Twilight Positioning
shown in Sandstorm



**Pocketed Coil
Spring Seat
Standard!**



TWILIGHT

**Orion
PR405-MLA**

Introducing the world's first three-position power recliner with Twilight Technology! Our dual-motor system comes standard and allows for a deeper recline to achieve rejuvenating positions no other three-position chair on the market offers. Rest, relax and rejuvenate in the Orion with Twilight for a weightless experience!

 **GOLDEN**

**PROUDLY MADE
BY AMERICANS**
with Domestic & Imported Materials


PR405-MLA

Model Number	PR405-MLA
Size	Medium/Large
Weight Capacity	375 LBS
Recline Positions	4-Way
Standard Fabrics	5
Fabric Options	All
Chaise Seat	YES
Coil Springs	YES
Back Type	Seam
Head Rest Cover	None
Arm Covers	None
Battery Back Up	YES
Overall Width	40.5"
Overall Height	46.5"
Width Between Arms	23"
Floor to Top of Seat	21"
Seat Depth	21"
Seat to Top of Back	30.5"
Seat to Top of Arm	5"
Seat Angle Raised	17°
Floor to Top of Foot Rest Reclined	25"
Floor to Top of Back Reclined	35.5"
Back Angle Closed	105°
Back Angle Reclined	145°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	21"
Extension reclined from Back of Seat to end of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Option Available	NO
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	YES
Extra Pocket Option Available	Standard
Heat & Massage Option Available	NO
Gel Memory Foam Option Available	NO
Power Pillow Option Available	NO
AutoFlex Option Available	NO
Foot Rest Extension Option Available	NO
Straight Lift Option Available	NO
USB Port	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.




TWILIGHT

**Motor #1**

Controls the movement of both the backrest and footrest.

Motor #2

Activates Twilight technology, a tilting movement that provides a weightless experience.

**Smart Tek**

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

**Battery Backup System**

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

