


M A X I C O M F O R T TM S E R I E S



**Now with
Power Headrest
& Power Lumbar
STANDARD!**



Hazelnut



Calypso



Shiraz



Coffee Bean*

* upcharge applies

Day Dreamer
PR632
shown in Coffee Bean



PR632

The DayDreamer™ is a medium sized power lift recliner offering our exclusive PowerPillow™ articulating headrest and power lumbar STANDARD! It features The Cloud™ seating system and our custom bucket seat design. Recline, relax, and adjust your head and back positions with the touch of a button.





Day Dreamer

Model Number	PR632
Size	Medium
Weight Capacity	375 Lb.
Positioning	MaxiComfort
Chaise Seat	YES
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Option	NO
Head Rest Cover	NO
Arm Rest Covers	NO
Battery Back Up	YES
Overall Width	38.5"
Overall Height	45"
Width Between Arms	22.5"
Floor to Top of Seat	20"
Seat Depth	20.5"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	24"
Distance Required from Wall Reclined	22"
Extension reclined from Back of Seat to End of Foot Rest	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option	YES
Extra Pocket Option	Standard
Heat & Massage Option	YES
Power Pillow Option	Standard
Power Adjustable Lumbar Option	Standard
AutoFlex Option	NO
Foot Rest Extension Option	NO
Straight Lift Option	NO
USB Charging Port	NO
Place of Origin	USA



Power Articulation Head Forward



Power Articulation Head



Power Articulation Head Back



Power Adjustable Lumbar

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

