

S I G N A T U R E S E R I E S



Space Saver
PR931 M



Anchor



Oxford



Bittersweet



Sandstorm



Port



Are you short on space, but searching for extreme comfort? You don't have to look any further than our space saving power lift recliner. With a seamed back and forward-gliding recline motion, this design offered in medium requires only 5" inches from any wall. The Golden Space Saver will fit you perfectly in almost any sized room!

PR931





PR931

S I G N A T U R E S E R I E S

	Space Saver
Model Number	PR931M
Size	Medium
Weight Capacity	350 LBS
Positioning	Space Saver
Chaise Seat	YES
Pocketed Coil Seat	NO
Back Type	Seam
Removable Back Option	YES
Head Rest Cover	YES
Arm Rest Covers	NO
Battery Back Up	YES
Overall Width	32.5"
Overall Height	41.5"
Width Between Arms	20"
Floor to Top of Seat	19"
Seat Depth	21"
Seat to Top of Back	28.5"
Seat to Top of Arm	5.5"
Seat Angle Raised	20°
Floor to Top of Foot Rest Reclined	22"
Distance Required from Wall Reclined	5"
Extension reclined from Back of Seat to End of Foot Rest	42.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option	NO
Extra Pocket Option	NO
Heat & Massage Option	YES
Power Pillow Option	NO
Power Adjustable Lumbar Option	NO
AutoFlex Option	NO
Foot Rest Extension Option	NO
Straight Lift Option	NO
USB Charging Port	YES
Place of Origin	USA



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

