

M I D S I Z E S C O O T E R



Seat Color  
Black Vinyl

Body Colors  
Red

Blue

# GOLDEN COMPANION

3 wheel (Midsize)  
GC240C

The Companion scooter provides great style with high performance and user-friendly disassembly. New features include an adjustable LED headlight; new ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion offers great value in a full-time portable scooter.



### Companion 3 wheel Specifications

<b>Model</b>	<b>GC240C</b>
Medicare Code	K0807
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5.3 mph
Operating Range <sup>1</sup>	19.9 miles
Ground Clearance Under Deck	3.75"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	37"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
<b>Dimensions:</b>	
Length	42.25"
Width	23.5"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	23" - 25"
Ground to Top of Deck	6.25"
Deck to Top of Seat	17.5" - 19.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
<b>Weight of Unit: (assembled)</b>	163 lbs
Front Half	45 lbs
Rear Half	48 lbs
Seat w/Arms	26 lbs
Batteries	21 lbs ea
<b>Tires:</b>	Solid, Low Profile
Front	9"
Rear	10.75"
Rear Anti-Tippers	3"
<b>Standard Seat:</b>	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" x 16"
Color	Black Vinyl
<b>Optional Seats:</b>	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black Vinyl
<b>Seat Options:</b>	
Pneumatic Rise Seat	Yes
Power Elevating Seat	Yes
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years
Batteries	1 Year



Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment



Patented infinite adjustable tiller




Quiet & maintenance free motor

### Standard Features on Every Model:

- Under mat troubleshooting guide
- Delta tiller with auto lock mechanism
- Easy reach tiller adjustment
- Adjustable sliding seat
- Adjustable LED headlight
- U-1 batteries standard



Disassembles easily!



<sup>1</sup> Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

<sup>2</sup> Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

